

Micro Centrifuges

MCF08VT / HMCf12VT

LabQuest micro centrifuges are engineered for advanced aerodynamics to ensure low heat generation due to air friction and low noise operation.

Safety features include imbalance detection with auto cut-off. Electronic brake for auto stop on lid opening. Rotor and adaptors are supplied as standard accessories.

**Accessories
FREE**





MCF08VT

Code : **100MCF00800VT0**



HMC12VT

Code : **100HMC1200VT0**

Technical Information

Specification	MCF08VT	HMC12VT
Rotor position	8	12
Tube capacity in ml	1.5 / 2.0 / 0.2 / 0.5	1.5 / 2.0 / 0.2 / 0.5
Max rpm	6000	15000
RCF	2000	15596
Display	digital	digital
Noise level	< 55db	< 60db
Imbalance detection	No	Yes
Dimension (W X D X H) in mm	162 x 157 x 116	260 x 244 x 205
Product code	100MCF00800VT0	100HMC1200VT0

BOROSIL®

Pioneers in speciality glass

www.borosil.com

Borosil Glass Works Ltd.

Corporate Office

1101, Crescenzo, G-Block, Opp. MCA Club, Bandra Kurla Complex, Bandra (E), Mumbai - 400051, India
Tel : +91 22 6740 6300 • Fax : +91 22 6740 6514 • Email : borosil@borosil.com

For export enquiries write to us eximbor@borosil.com