

**Manufacturer, Exporter & Supplier of Scientific Instruments,
Microscope, Centrifugal Machines, Pharmacy Equipments, Clean Air Equipments,
Physical Lab & Analytical Instruments**



STUDENT MICROSCOPE

Model : LS-801 STUDENT MICROSCOPE

MODEL NO.	PRICE RS.
Model : LS-801	

SPECIFICATION :

Body : Sturdy and well balance aluminum casted body inclinable to 90 Degree.

Stage : Fix square stage 110 x 110mm.



Nose Piece : Triple Revolving nose piece with smooth click stops.

Focusing : Focusing by coarse & fine focusing knob.

Condenser : Sub stage abbe type N.A. 1.25 condenser with IRIS Diaphragm (Fix).

Optics : Parafoled Achromat objectives 10 x 45 x(SL) & Eye Piece 10 x Huygenian.

Illumination : Illumination through the reflecting mirror.





STUDENT MICROSCOPE




Model : LS-801B STUDENT MICROSCOPE (All Brass)

SPECIFICATION :

- Body** : Sturdy and well balance aluminum casted body inclinable to 90 Degree.
- Stage** : Fix square stage 110 x 110mm.
- Nose Piece** : Triple Revolving nose piece with smooth click stops.
- Focusing** : Focusing by all Brass coarse & fine focusing knob.
- Condenser** : Sub stage abbe type N.A. 1.25 condenser with IRIS Diaphragm (Fix).
- Optics** : Parafoiled Fully Brass Achromat objectives 10 x 45 x(SL) & Eye Piece 10 x Huygenian.
- Illumination** : Illumination through the reflecting mirror.

MODEL NO.	PRICE RS.
Model : LS-801B	

Optional: LED Lamp



**ALL BRASS
HEAVY BASE**