



Product Catalogue

MICROSCOPES & ACCESSORIES

#magnusmicroscopes



ABOUT US

Driven and committed to developing and manufacturing products for today's scientific and medical community. We strive to continuously evolve to stay abreast with the speed of change in the technology world.

CH20i

Biological Microscope

OLYMPUS

- Anti-fungus optics
- Optics with multi-layer coating
- Easily Interchangeable objectives
- Universal power supply 100V ~ 240V
- Main operating controls within easy reach
- Optimization of illumination with aspheric lenses
- Available in monocular, binocular and trinocular version



Model	Head	Objectives	Eyepiece	Stage	Condenser	Illumination
CH20i Binocular	Binocular (inclination 45°)	iNEA achromat 4x, 10x, 40x, 100x	iCWHK 10x (FN. 18)	Mechanical stage 120mm X 132mm	1.25 Abbe (oil immersion)	6v 20w Halogen
CH20i Trinocular	Trinocular (inclination 45°)					

