

H



- ✦ Measuring Range - $0.000 \sim 3.0000 \text{ g/cm}^3$
- ✦ Data Storage - 32GB
- ✦ Resolution - 0.001 g/cm^3
- ✦ Calibration - Dry Air / Distilled Water
- ✦ Temperature Control Range - $15^\circ\text{C} \sim 40^\circ\text{C}$
- ✦ Resolution - 0.01°C
- ✦ Observation - High Definition Camera / Video Observation, Density g/cm^3 , Density kg/cm^3 (Specific Gravity of Water)
- ✦ Measurement Mode - Fast Pre-deduct / Accurate Measure

Automatic Density Meter

ADM3